

PLAIN DRAPE



DESCRIPTION

Standard drape sheet, reinforced

Item#	size
D110	155 x 250 cm

TOP/SIDE/BOTTOM DRAPE



DESCRIPTION

Absorbent reinforcement, (74 x 25 cm)
Hook-and-loop line holders

Item#	size
D120	155 x 250 cm

TABLE COVER




DESCRIPTION

Standard back table cover,
SMS reinforcement

Item#	size
D131	155 x 130 cm
D132	155 x 150 cm
D133	155 x 200 cm
D134	155 x 250 cm

MAYO STAND COVER




DESCRIPTION

Mayo stand cover, reinforced

Item#	size
D140	78 x 145 cm

ANGIOGRAPHY DRAPE




DESCRIPTION

Absorbent, impervious reinforcement 60x75cm
Fenestration, dual circular

Item#	size
D150	155 x 300 cm

ANGIOGRAPHY WITH POUCH




DESCRIPTION

Absorbent, impervious reinforcement 60x75cm, Pouch 30X50cm
Fenestration, dual circular

Item#	size
D160	155 x 300 cm

SPLIT U-DRAPE




DESCRIPTION

- Absorbent reinforcement (60 x 75 cm)
- Split with adhesive (15 x 45cm)
- Hook-and-loop line holder

Item#	size
D171	Leg-U (155x 210cm)
D172	Arm-U (155x 175cm)

ANCILLARY DRAPE




DESCRIPTION

Ancillary drape sheet, (60x75cm) reinforced fenestration

Item#	size
D181	155x 110 cm
D182	155x 300 cm

CRANIOTOMY DRAPE



DESCRIPTION

Absorbent, impervious reinforcement, 60x 75 cm.
Fenestration, oval with full incise
Fluid control pouch 50x 50 cm


Item#

size

D190

300x 320 cm

NEUROLOGICAL DRAPE



DESCRIPTION

Absorbent, impervious reinforcement, 60x 75 cm
Fenestration, oval with full incise
Fluid control pouch 50 x50 cm


Item#

size

D200

155 x 155 cm

ARTHOSCOPY DRAPE



DESCRIPTION

Absorbent, impervious reinforcement, 60x 75 cm.
Fluid control pouch 50x 50 cm

Item#

size

D210

155x 300cm

SPINE DRAPE



DESCRIPTION

Absorbent, impervious reinforcement, 60x 75 cm
Fenestration, rectangle with full incise


Item#

size

D220

155 x 300 cm

LAPAROTOMY DRAPE



DESCRIPTION

Absorbent, impervious reinforcement, 60 x 75 cm.
Cord-holding tabs

Item#

size

D230

155 x 280 cm

CARDIO VASCULAR DRAPE



DESCRIPTION

Absorbent, impervious reinforcement
Split with adhesive

Item#

size

D240

300x 320 cm

CYSTOSCOPY T-DRAPE



DESCRIPTION

Absorbent reinforcement, 60x 75 cm.
Cord-holding tabs
Fenestration


Item#

size

D250

300x 320 cm

EXTREMITY DRAPE



DESCRIPTION

Absorbent, impervious reinforcement, 60x 75 cm.
Fenestration, oval with full incise

Item#

size

D260

155 x300 cm

CUT SHEET



DESCRIPTION

Blue Reinforced


Item#

D270

size

75 x200 cm

EYE DRAPE



DESCRIPTION

Eye drape with fluid collection pouch

Item#

D280

size

75x 75 cm

TROLLEY COVER



DESCRIPTION

Poly Vinyl


Item#

D290

size

130x 175 cm

DRAPE TOWEL & WHITE TOWEL



DESCRIPTION

Blue drape towel, absorbent
with adhesive strip, sterile
White absorbent towel

Item#

D301

D302

D303

size

40x 75 cm

80x 75 cm

White (27x 28cm)